



Product description

Order number	55.495
Product description	Tube, 8.5 ml, 75 x 15.7 mm, round base, PS, 1,000 pcs./bag

Product characteristics

Graduation	No
Type of cap	Cap
Base shape	Round base
Centrifugation max (RCF)	4000 x g

Size

Sample volume	8,5 ml
Diameter	15,7 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polystyrene (PS)
Colour of product	Transparent
Tube material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	Sealed at 0.6 bar/10 min.
Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

Packaging

Minimum order qty.	2000
Type of packaging	Bag
Pcs. / inner box	1000
Pcs. / box	2000
Pcs. / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0,0539 cbm
Weight of case	9,08 KG
EAN of inner box	4038917082013
EAN of outer box	4038917001120

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

