



Product description

Order number	32.329
Product description	Tube 4.5 ml, Clotting Activator/Serum, 75x13 mm, white EU code, PP, print

Product characteristics

Type of preparation	Clotting activator
Colour code	EU
Label/ Print	With printed label
Colour of print/label	White
Graduation	Yes
Cap	Push cap
Type of cap	Cap
Base shape	Round base
Recommended centrifugation	2000 x g
Centrifugation time	10 min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	4,5 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	Low Density Polyethylene (LD-PE)
Colour of cap	White
Tube material	Polypropylene (PP)

Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	batched

Packaging

Minimum order qty.	1000
Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	64000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0,0255 cbm
Weight of case	4,23 KG
EAN of inner box	4038917081160
EAN of outer box	4038917000307

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

